



Neutral pH Non-Enzymatic Low Foam Detergent Clean Start Chemicals

Specifically designed for Skytron washer-disinfector systems.

Scan the QR code
to see Skytron's full
offering of chemicals.



Neutral pH Non-Enzymatic Low Foam Detergent

Cleaning and Disinfection Solution

Neutral pH Non-Enzymatic Low Foam Detergent, a highly concentrated neutral detergent, is recommended for use on all types of instruments and equipment, particularly for those that have already been pre-cleaned. With built-in lubrication agents, it not only ensures safety but also prolongs the life of instruments.

Recommended Uses:

- Medical/surgical instruments
- Stainless steel and aluminum items
- Respiratory therapy equipment
- Evacuation systems
- Tungsten and carbide copper items

Technical Data:

- pH neutral
- Non-corrosive
- No phosphates, chlorides or hydroxides
- Dilution – 1/8 to 1 oz. per gallon of water

Advantages:

- Safe and effective for stainless steel, aluminum, tungsten, carbide copper, and brass
- Improves cleaning and rinses freely
- Formulated for all brands of automatic washers and sterilizers
- Built-in lubrication agents to prolong instrument life

Packaging:

- Four 1 gallon containers – CSTART-9601
- Two 2.5 gallon containers – CSTART-9625
- One 5 gallon container – CSTART-9605

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV2 08-24
PN L2-010-250