



Mobile Surgical Supply Cabinets

Skytron Stainless

Creating the perfect clinical environment starts with top-notch storage and organization solutions. Skytron offers seamless functionality for your healthcare setting. Made in the USA with unmatched quality to reduce upkeep, our cabinets stand out for their cost efficiency and long-lasting performance.



Mobile Surgical Supply Cabinets



Standard Features

- Framed double-pan glazed doors hung on full length stainless steel hinges
- 1/4" tempered safety glass
- Top and bottom spring loaded nylon roller door catches
- Stainless steel shelves adjustable in 1/2" increments
- Fixed bottom shelf
- Copper pull handles

Optional Features

- Additional shelves
- Pull-out cath hook system to hang inventory
- Scope hangers for easy and safe endoscope storage
- Adjustable height hanger hooks



Scope hangers for endoscope storage



Pull-out cath hook system

Single Door Cabinets	
Model	MSSC246024
Exterior Dimensions	
Width (A):	26" (661 mm)
Depth (B):	28" (711 mm)
Height (C):	68" (1727 mm)
Interior Dimensions	
Width (D):	21" (533 mm)
Depth (E):	23" (584 mm)
Height (F):	57" (1448 mm)
Weight	
Empty:	201 lbs (91 kg)
Loaded:	401 lbs (182 kg)
Crated:	246 lbs (112 kg)

Double Door Cabinets				
Model	MSSC306024	MSSC366024	MSSC426064	MSSC486024
Exterior Dimensions				
Width (A):	32" (813 mm)	38" (965 mm)	44" (1118 mm)	50" (1270 mm)
Depth (B):	28" (711 mm)	28" (711 mm)	28" (711 mm)	28" (711 mm)
Height (C):	68" (1727 mm)	68" (1727 mm)	68" (1727 mm)	68" (1727 mm)
Interior Dimensions				
Width (D):	27" (686 mm)	33" (838 mm)	39" (991 mm)	45" (1143 mm)
Depth (E):	23" (584 mm)	23" (584 mm)	23" (584 mm)	23" (584 mm)
Height (F):	57" (1448 mm)	57" (1448 mm)	57" (1448 mm)	57" (1448 mm)
Weight				
Empty:	210 lbs (95 kg)	247 lbs (112 kg)	282 lbs (128 kg)	318 lbs (144 kg)
Loaded:	410 lbs (186 kg)	445 lbs (202 kg)	512 lbs (232 kg)	518 lbs (235 kg)
Crated:	256 lbs (116 kg)	291 lbs (132 kg)	326 lbs (148 kg)	364 lbs (165 kg)

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV1 04-24
PN L2-010-252